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EN 61010-1
 EN 61010-2-051

Thank you for purchasing a Heidolph Instruments product. This item has been designed, made and inspected in compliance with DIN EN ISO 61010 for long-term and flawless operation.

SUMMARY

Summary	16
standard hardware & accessories	17
general information	18
Safety information	18
set-up	19
operation and controls	20
Turn item ON.....	20
Speed setting.....	20
Torque display.....	21
Interfaces.....	22
Protective equipment.....	24
cleaning & servicing	25
disassemble & storage	25
disposal	25
troubleshooting	25
specifications	26
warranty, liability & copyright	28
questions / repair work	28
CE-declaration of conformity	28

Important information

Advice about power cord / mains supply

Caution: mandatory action

Caution: fire- and explosion hazard

Advice about maintenance / repair

STANDARD HARDWARE & ACCESSORIES

Standard Hardware

Product	quantity	P/N	
		230 V / 50/60 Hz	115 V / 50/60 Hz
RZR 2051 control or	1	501-20511-00	501-20511-01
RZR 2052 control or	1	501-20521-00	501-20521-01
RZR 2102 control or	1	501-21021-00	501-21021-01
RZR 2102 control Z	1	501-21024-00	501-21024-01
support rod	1	22-02-14-01-08	22-02-14-01-08
2,5 mm Allen key	1	02-07-01-01-04	02-07-01-01-04
3,0 mm Allen key	1	02-07-01-01-02	02-07-01-01-02
10.5 mm chuck	1	509-02000-00	509-02000-00
10.5 mm chuck (only RZR 2102 control Z)	1	11-001-001-56-0	11-001-001-56-0
T-handle (for chuck)	1	02-07-03-01-02	02-07-03-01-02
Instruction Manual	1	01-005-002-95	01-005-002-95



Accessories & Stirring tools (optional)

Product	P/N
Stand	570-12000-00
Clamp	570-22000-00
protective sheath	509-08000-00
flexible stirrer coupling shaft	509-03000-00
interface extension (IN)	591-41000-00
remote control	591-40000-00
RS232- St15/Bu9 interface cable(stirrer – PC)	14-007-045-17
SUB-HDD 15 interface connector ass'y. (for connecting analog interface)	14-010-009-24
stirring paddle BR 10	509-10000-00
stirring paddle BR 11	509-11000-00
stirring paddle BR 12	509-12000-00
stirring paddle BR 13	509-13000-00
stirring paddle BR 14	509-14000-00
prop-type stirring paddle PR 30	509-30000-00
prop-type stirring paddle PR 31	509-31000-00
prop-type stirring paddle PR 32	509-32000-00
prop-type stirring paddle PR 33	509-33000-00
turbine TR 20	509-20000-00
turbine TR 21	509-21000-00
half moon HR 18	509-18000-00
anchor AR 19	509-19000-00
prop-type stirring paddle PR 39	509-39000-00
enclosure NS 29/32	509-09000-00



GENERAL INFORMATION

-  **Unpack your item carefully.**
Inspect for damage and report such damage or missing parts to your supplier right away.
-  Read your Instruction Manual carefully. Take time to save time when working with your product. Make sure that every user has read and understood the Instruction Manual.
-  Please store the Instruction Manual in a place easily accessible to every user.

IF ALL ELSE FAILS, READ THESE INSTRUCTIONS !

-  A so-called EURO-plug (DIN 49441 CEE 7/VII 10/ 16 A 250 V) is standard on all of the products. For the Continental US they feature a US-standard plug (NEMA Pub.No.WDI.1961 ASA C 73.1 . 1961 page 8 15A 125V).
-  For using the item in a country with deviating outlet / plug systems, we recommend to use approved adapters or to have an electrician replace the standard plug with one suiting your needs.
-  As shipped, the item features a protective ground wire. When replacing the original plug, make sure to reconnect this protective ground wire in the new plug !


SAFETY INFORMATION

-  **Please comply with all safety and accident-prevention regulations as in force for laboratory work !**
-  **Use extra care when working with flammable substances; refer to safety data sheets.**
-  **When connecting your item with your local power supply, please make sure your item is designed for your local voltage; go by the data plate on the item.**
-  **Turn your power switch OFF whenever the item is not used, or before disconnecting the plug.**
-  **Repair work is limited to technicians approved by Heidolph Instruments.**
-  **Use extra care when working in the vicinity of flammable and explosive substances. Motors are of non-sparking type, the item itself however is not explosion-proof.**
-  **Please connect your stirrer with a protective-ground outlet only.**
-  **Please strictly comply with current safety regulations for machines exposing rotating shafts.**
-  **Your item needs a solid stand.**

SET-UP

1. Install stand

make sure the indent (A) of the rod is in line with threaded pin (B). Install stand in hole (C) in the rear of the item. Secure (pull tight) threaded pin (B) using Allen key (D).



2. Install chuck

make sure the indent (E) in the drive shaft is in line with threaded pin (F) in the chuck. Secure threaded pin (F) with Allen key (H) in your hardware bag.

3. Install stirrer on stand

Install clamp (I) on stand (J).
Now introduce stirrer attaching rod (K) into clamp (I).
Adjust position and secure (pull tight) with two clamp screws (L).

4. Install stirring tool


Introduce stirring tool (N) in chuck (G), put in place and secure with T-handle (M). Immersion depth of your stirring tool (N) is adjusted by displacing the tool in and out in the open (G) chuck (except RZR 2102 control Z).
Refer to arrow mark.

5. Install flexible coupling shaft (P/N 509-03000-00)

Hold the flexible coupling shaft (O) in the chuck (G) by its metallic stud (P) and secure with T-handle. At this point introduce stirring tool's (N) shaft end in collet (Q) and pull tight (CW turn).

When using the flexible coupling shaft make sure to install an adequate run-out protection, e.g. our enclosure (P/N 509-09000-00)


6. Install shaft guard (P/N 509-08000-00)

A shaft guard (R) is available as an extra safety protection feature. This guard tube is held by 2 ea. screws (S) located at the top of the stirrer chuck.
This guard, attached to the stirrer housing, is used to protect the user from being hurt by the fast-spinning shaft.

OPERATION AND CONTROLS

1. Turn item ON

Set speed control (4) to its extreme left position (CCW). When turning the power switch (1) ON, the following info will appear in the display (5)

Heidolph
RZR V2.X

followed by speed and torque; at standstill, „rpm“ in line 1 will change between rpm and OFF, which means the item is ready to go.

2. Speed setting

2.1 RZR 2051 control, RZR 2052 control, RZR 2102 control and RZR 2102 control Z

At standstill of the item, speed may be preselected with the rotary speed control (4). RZR 2051 control speed range includes 50 – 2000 rpm, the RZR 2051 control features 30 – 1000, the RZR 2102 control features 12 – 400 and 60 – 2000 in 2 ranges and the RZR 2102 control Z features 4 – 108 and 17 – 540 in 2 ranges. Speed selected appears in line 1 of display (5).

Hitting the ON/OFF key (2), the stirrer is going to speed-up to the level selected, after a short pause.

Acceleration is smooth, thus preventing the medium from splashing.

Whenever changing speed with rotary control (4) while running, display is going to show nominal speed.

For convenience, nominal speed is displayed with the "⊖" prefix. This prefix disappears, when actual speed is displayed !

2.2 RZR 2102 control and RZR 2102 control Z

This stirrer features two different speed ranges:

RZR 2102 control:

range I: 12 – 400 rpm

range II: 60 – 2000 rpm

RZR 2102 control Z:

range I: 4 – 108 rpm

range II: 17 – 540 rpm

Transmission ratio is obtained with the **RANGE** key (8). Speed range actually selected is displayed by an LED (9 or 10).

At standstill, speed may be selected with the speed control (4), within both speed ranges (I or II). In standby mode, the display (5) shows the speed value (set-mode).

Change speed range, RZR 2102 control and RZR 2102 control Z


Select your speed range with the **RANGE** key (8). LEDs (9) and (10) confirm the range selected. When switching from range I to II, stirring speed is multiplied by factor five.

-  **Note:**
To avoid splashes (in particular when open containers are used), speed range cannot be changed unless the item has been turned OFF before.

How to proceed:

- turn stirrer OFF with ON/OFF key (2)
- turn speed control (4) to its extreme left position (CCW) (= low speed)
- change speed with key (8)
- turn item ON with ON/OFF key (2)
- set speed with speed control (4)
- speed appears in display (5), line 1

You may change from range II to range I at any time without any extra precautions necessary.

3. Torque display

Both stirrers feature a digital display (5) for stirring torque (line 2). Before measuring, you need to „zero“ your torque display (= 0,0 Ncm), to benefit of the highest accuracy available from this system (refer to stipulations below). This is of particular importance, when changing stirrer speed.

How to „zero“ the torque display:

- allow the stirrer to warm at max. speed for about 20 minutes.
- define speed needed for your particular job
- select standby mode with the ON/OFF key (2)
- lift stirring tool with stirrer out of the medium, or remove stirring tool from chuck (refer to para. 4 „set-up“). This way, the motor won't be able to build up any torque while spinning.
- turn item ON with ON/OFF key (2)
- with the motor running, hit the „cal“ key (7), display (5) is going to change to „0,0“ Ncm in line 2
- turn the item OFF, re-install stirring tool in (refer to para. 4 „set-up“), or immerse your system in the medium again
- at this point you have an accurate torque reading

4. Interfaces

4.1 Interface connection

Turn power switch (1) to OFF before attempting to connect an interface.

To connect your interface, you will need an interface cable with SUB-HDD 15 connector which fits the socket (11) located in the item's rear panel.

For more information about wiring order, refer to „analog interface and serial RS 232 interface“ diagram.

Only connect but the wires shown to avoid malfunction of your interface

Serial RS 232 interface

stirrer end


PC end (COM x)

SUB-HDD 15 plug

e.g. Mc Murdo HDE15PTD

SUB-D 9 socket

e.g. Mc Murdo SDE9ST

P/N (options): 14-007-045-17-0 RS232- St15/Bu9 interface cable

analog interface 0 – 10 V


stirrer end

user-defined end

SUB-HDD 15 plug

e.g. Mc Murdo HDE15PT

plug connectors

P/N (options): 14-010-009-24-0 SUB-HDD 15 plug ass'y.

4.2 Serial RS 232 interface

This type of interface allows to control your stirrer with a PC. For connecting the serial interface, please refer to 4.1.

Interface commands

OUT		IN	
command	action	command	action
r	speed [rpm] OUT –	Rxxxx	select speed (1-4 digits) (enter with CR-key)
s	nominal speed [rpm] OUT	A	range II
m	Torque [Nmm] OUT	B	range I
f	display error status: No Error! Motor Error! Motor Temperature!	C	delete error mode- motor restart
		D	RC OFF, motor control uses rheostat
		N	set „0“-reference mark for speed OUT information

Interface parameter: 9600 baud, no parity, 8 bit, 1 stop bit

4.3 Analog interface 0 – 10 volts

OUT

The analog 0-10 volts interface sends torque and speed information. For connecting the analog interface, please refer to 4.1.

OUT:	RZR 2051	0 rpm = 0V 50 – 2000 rpm = 0,25V – 10V	0 – 20 Ncm = 0V – 10V
	RZR 2052	0 rpm = 0V 30 – 1000 rpm = 0,3V – 10V	0 – 100 Ncm = 0V – 10V
	RZR 2102	0 rpm = 0V	
	<i>Range I</i>	12 – 400 rpm = 0,3V – 10V	0 – 200 Ncm = 0V – 10V
	<i>Range II</i>	60 – 2000 rpm = 0,3V – 10V	0 – 40 Ncm = 0V – 10V
	RZR 2102 Z	0 rpm = 0V	
	<i>Range I</i>	4 – 108 rpm = 0,3V – 10V	0 – 740 Ncm = 0V – 10V
	<i>Range II</i>	17 – 540 rpm = 0,3V – 10V	0 – 150 Ncm = 0V – 10V

Interface extensions (IN); P/N 591-41000-00

refer to options section

Using an interface extension card, nominal speed of the stirrer can be triggered at 4 – 20 mA.

RC (remote control) (P/N 591-40000-00)

An ON/OFF remote control is optional.

Remote control allows ON/OFF mode of the stirrer (same as On/Off key)

5. Protective equipment

5.1 Overload LED

The stirrer features and extra overload LED (6) that will light in an overload situation. To get higher torque out of the item, overload operation is acceptable for a short period of time. After extended overload operation, a built-in thermal circuit breaker will disconnect the item automatically.

the overload LED (6) must not light during continuous operation; the item will be cut-off automatically after a short time!

5.2. Thermal overload protection

The stirrer features built-in overheat protection. If the item is getting too hot, either by overload (refer to 5.1) or high ambient temperatures, it will be disconnected automatically. The following error prompt appears in the display (5):

XXX rpm/off
OV-Temp!

Allow item to cool-off. As soon as temperature drops to an acceptable level, display (5) will read:

XXX rpm/off
Mot-Err!

Hitting the ON/OFF key (2), the item will return to normal operation, running at speed set.

5.3. Mechanical overload protection

In case of an overload situation including excessive torque (refer to torque graph under "Specifications"), or by stuck stirrer shaft, the stirrer is going to disconnect automatically.


Display (5) reads:

XXX rpm/off
OV-Temp!

After 90 seconds, display (5) reads:

XXX rpm/off
Mot-ERR!

When having fixed the problem, hit ON/OFF key (2) for re-start. Item will return to normal operation, running at speed set.

CLEANING & SERVICING

Cleaning: wipe housing clean with a damp cloth (add some sort of mild liquid soap).

Note

To avoid damage to the surface finish, avoid using chlorine bleach, chlorine-based detergents, abrasive substances, ammonia, rags or cleaning agents containing metal particles.

The item is maintenance-free. Repair work is limited to technicians approved or appointed by Heidolph Instruments. Please call your local Heidolph Instruments Dealer or a Heidolph Instruments Field Representative (also refer to page 28)

DISASSEMBLE & STORAGE

Disassemble (refer to page 19, „set-up“)

Turn the item OFF and pull power plug.

To disassemble the item easily, remove all flasks, containers or other appliances in the vicinity of the stirrer.

Open chuck with T-handle and remove stirring tool.

Unlock clamp and remove stirrer head from head clamp.

With 2,5 mm Allen key loosen threaded pin and remove chuck. With the 3 mm Allen key loosen and remove stand rod.

Forward & Store

1. We recommend to store the item and its components in its original box, or a similar container that offers adequate protection against damage in transit. Tape the box securely.
2. Store the item in a dry place.

Caution

Do not jolt or shake the item during transport.

DISPOSAL

For disposal, please comply with your local or national regulations.

Split by metal, plastic, etc.

Packing material to be treated as described above (material split).

TROUBLESHOOTING





1. power switch dim
 - check power supply and power cord
2. won't stir:
 - safety installation intervened, refer to display reading (refer to 5 „safety installations“)
 - interface connection miswired (refer to 4 „interfaces“)
3. speed selector (range II) won't work:
 - will work at standstill only (refer to 2 „speed setting“)

In case of a malfunction and you cannot correct with above instructions, contact your authorized Heidolph Instruments Dealer.

SPECIFICATIONS

power dissipation (nom. load) RZR 2051 control	P = 80 watts
power dissipation (nom. load) RZR 2052 control	P = 140 watts
RZR 2102 control	P = 140 watts
RZR 2102 control Z	P = 140 watts
power rating, motor, RZR 2051 control	P = 50 watts
power rating, motor, RZR 2522 control	P = 100 watts
RZR 2102 control	P = 100 watts
RZR 2102 control Z	P = 100 watts
dimensions (in mm) RZR 2051 control	72 wide; 206 high; 176 deep
dimensions (in mm) RZR 2052 control	82 wide; 211 high; 176 deep
dimensions (in mm) RZR 2102 control Z	82 wide; 292 high; 176 deep
weight RZR 2051 control	2.8 kg
weight RZR 2052 control	3.7 kg
weight RZR 2102 control	3.7 kg
weight RZR 2102 control Z	4.7 kg
operating temperature	0 – 40 °C
storage temperature	- 20 to + 80 °C
relative humidity	95 % (avoid reaching dewpoint)
protection class	IP 40
Torque	refer to graphs below
speed range	
RZR 2051 control	50 – 2000 rpm
RZR 2052 control	30 – 1000 rpm
RZR 2102 control	step 1: 12 – 400 rpm step 2: 60 – 2000 rpm
RZR 2102 control Z	step 1: 4 – 108 rpm step 2: 17 – 540 rpm
serial interface	RS 232
analog interface, (OUT), 0 – 10 volts	torque and speed
interface extension (IN) 4 – 20 mA; 0 – 10 volts	speed and speed range (optional)
remote control	ON/OFF (optional)
stirrer coupling RZR 2051 control, RZR 2052 control and RZR 2102 control	10,5 mm chuck, drive shaft features bore to hold tool shafts of up to 10 mm dia.
stirrer coupling RZR 2102 control Z	10,5 mm chuck, drive shaft don't features bore to hold tool shafts.
stirrer attaching system	13 mm dia. Rod (stand)
motor protector	thermal circuit breaker disconnects motor in an overload situation (permanent cut-off); RESET with ON/OFF key
pilot light, „item ON“	power switch (GREEN)
shaft guard	fits shaft guard (optional)
BCI/TVI suppression	complies with EN 50081
operational mode	100% continuous
display	8-digit LCD, 2 lines, 5.6 mm high characters
measuring range, display	4 – 2000 rpm / 0-700 Ncm
Resolution	1 rpm / 0,1 Ncm

Torque graphs

WARRANTY, LIABILITY & COPYRIGHT

Warranty

Heidolph Instruments guarantees warranty that the present product shall be free from defects in material (except wear parts) and workmanship for 3 years from the date shipped off the manufacturer's warehouse.

Transit damage is excluded from this warranty.

To obtain such warranty service, contact Heidolph Instruments (Tel.: +49 - 9122 - 9920-68) or your local Heidolph Instruments Dealer. If defects in material or workmanship are found, your item will be repaired or replaced at no charge.

Misuse, abuse, neglect or improper installation are not covered by this warranty guarantee..

Alterations to the present warranty guarantee need Heidolph Instruments' consent in writing.

Exclusion Clause

Heidolph Instruments cannot be held liable for damage from improper use or misuse. Remedy for consequential damage is excluded.

Copyright

Copyright in pictures and wording of the present Instruction Manual is held by Heidolph Instruments.

QUESTIONS / REPAIR WORK

If any **aspect** of installation, operation or maintenance remains unanswered in the present Manual, please contact the following address:

For repair services please call Heidolph Instruments (phone: +49 - 9122 - 9920-68) or your local, authorized Heidolph Instruments Dealer.

Note

You will receive approval for sending your defective item to the following address:

Heidolph Instruments GmbH & Co. KG
Lab Equipment Sales
Walpersdorfer Str. 12
D-91126 Schwabach / Germany
Tel.: +49 – 9122 - 9920-68
Fax: +49 – 9122 - 9920-65
E-Mail: sales@heidolph.de

Note

If you are based in the United States of America, please contact Heidolph US:

Heidolph Instruments LLC
Lab Equipment Sales
2615 River Rd.
Cinnaminson, NJ 08077
Phone: 856-829-6160
Fax: 856-829-7639
E-Mail: heidolph@snip.net

Safety Information

When shipping items for repair that may have been contaminated by hazardous substances, please:

- advise exact *substance*
- take proper *protective* measure to ensure the safety of our receiving and service personnel
- *mark* the pack IAW Hazardous Materials Act

CE-DECLARATION OF CONFORMITY

We herewith declare that the present product complies with the following standards and harmonized documents:

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© Heidolph Instruments GmbH & Co. KG