LABOROTA 4000eco

Betriebsanleitung
Instruction Manual
Mode d'Emploi
Instrucciones de
Empleo
Instruzioni per l'uso

D	DEUTSCH	Seite	3 - 35
E	ENGLISH	page	36 - 67
F	FRANCAISE	page	68 - 100
ES	ESPAÑOL	página	101 - 132
I	ITALIANO	pagina	134 - 167

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EN 61000-4-11:1994 + A1:2001

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EN 61010-1 + EN 61010-2-010

Thank you for purchasing a Heidolph Instruments product. This item has been designed, made and inspected in compliance with DIN EN ISO 61010 for long performance and continuous operation.

SUMMARY

STANDARD ITEM PLUS OPTIONS	38
GENERAL	39
SAFETY INFORMATION	39
SET-UP.....	41
1 Item set-up	41
2 Electric wiring	41
3 Install condenser support bar, when using G3 and G5eco glassware kits	43
4 Install vapor tube	44
5 Install glassware set G1	45
6 Install glassware set G3	47
7 Install glassware set G5eco	48
8 Cooling water fitting on glassware kit (G1, G3)	50
9 Connect vacuum pressure to glassware kit.....	50
10 Install Woulff flask / vacuum valve	51
11 Vacuum lines	52
12 Install protective hood, (option P/N 591.00010.00)	53
13 Install protective shield, (option P/N 591.00020.00)	54
CONTROLS AND OPERATION.....	54
14 Adjust evaporator flask angle	54
15 Use of hand lift	55
16 Adjust immersion depth of evaporator flask	55
17 Evaporator flask draw-off fixture.....	57
18 Displace heating bath.....	57
19 Fill hot bath	58
20 Work with protective hood and protective shield (optional, P/N 591.00010.00 and 591.00020.00)	59
21 Add substance to be distilled	59
22 Aerate by hand.....	59
23 Adjust hot bath temperature.....	59
24 Turn on and set rotary speed	60
INFORMATION ABOUT DISTILLATION SETTINGS	61
25 Rotary speed of evaporator flask	61
26 Temperature difference between hot bath and vapor.....	61
27 Temperature difference between vapor temperature and coolant.....	61
28 Solvent data	61
28.1 At boil temperature of 40°C.....	61
28.2 At boil temperatures deviating from 40°C.....	61

28.3 Solvents not depicted.....61

CLEANING AND MAINTENANCE.....63

DISASSEMBLE, FORWARD & STORE63

DISPOSAL64

TROUBLESHOOTING64

SPECIFICATIONS.....65

WARRANTY, LIABILITY & COPYRIGHT65

FAQ / REPAIR WORK66

CE-DECLARATION OF CONFORMITY.....67

Important information

Advice about power cord / mains supply

Caution: mandatory action

Caution: fire- and explosion hazard

Caution: hazard of burns, hot surface

Advice about maintenance / repair

STANDARD ITEM PLUS OPTIONS

	item	qty.	P/N 230/240V 50/60Hz	P/N 115V 50/60Hz
	basic item LABOROTA 4000eco	1	519-00000-00	519-00000-01
	heating bath WB eco	1	517-01000-00	517-01000-01
	glass set G1	1	513-00100-00	513-00100-00
or	glass set G3	1	513-00300-00	513-00300-00
	glass set G5eco	1	513-00501-00	513-00541-00
or	glass set G1 coated	1	513-00110-00	513-00110-00
or	glass set G3 coated	1	513-00310-00	513-00310-00
or	glass set G5eco coated	1	513-00511-00	513-00551-00
	loose parts:			
	coupling	1	23-09-03-01-03	23-09-03-01-03
	cocking spring	1	22-03-02-01-05	22-03-02-01-05
	PTFE 26 seal	1	23-30-01-01-30	23-30-01-01-30
	clamping sleeve 26	1	23-30-01-02-89	23-30-01-02-89
	coupling ring, vapor tube	1	11-300-005-19	11-300-005-19
	power cord	1	14-007-003-81	14-007-045-10
	link wire to heating bath	1	14-007-045-08	14-007-045-11
	Owner's Manual	1	01-005-002-73	01-005-002-73




Accessories

	item	P/N 230/240V 50/60Hz	P/N 115V 50/60Hz
	protective hood	591-00010-00	
	protective shield	591-00020-00	
	extension panels	591-00090-00	
	ROTAVAC valve control vacuum pump	591-00130-00	591-00130-01
	VAC control automatic vacuum controller	591-00340-00	591-00340-01
	vacuum switchbox for 3 ea. separate items	591-00400-00	591-00400-01
	ROTACOOOL	591-00200-00	591-00200-01
	Woulff-flask	591-00071-00	
	vacuum valve	591-24000-00	
	condensate cooler for ROTAVAC valve control	591-00081-00	



GENERAL

-  Unpack your item carefully. Inspect for damage and report such damage or missing parts to your supplier immediately.
-  Read your Instruction Manual carefully. Take the time to read the manual which will save time while working with your product. Make sure that every user has read and understood the Instruction Manual.
-  Please store the Instruction Manual in a place easily accessible to every user.

IF ALL ELSE FAILS, READ THESE INSTRUCTIONS !

-  A so-called EURO-plug (DIN 49441 CEE 7/VII 10/ 16 A 250 V) is standard on all of the products. For the Continental US they feature a US-standard plug (NEMA Pub.No.WDI.1961 ASA C 73.1 . 1961 page 8 15A 125V).
-  For using the item in a country with deviating outlet / plug systems, we recommend to use approved adapters or to have an electrician replace the standard plug with one mating your local system.
-  As shipped, the item features a protective ground wire. When replacing the original plug, make sure to reconnect this protective ground wire in the new plug !

SAFETY INFORMATION

-  **Please comply with all safety and accident prevention regulations as in force for laboratory work !**
-  **Use extra care while working with rotary evaporators.**
-  **Use eye protection and adequate clothing.**
-  **Use extra care while working with flammable substances; refer to safety data sheets.**
-  **Use extra care while working close to flammable and explosive substances. Motors are of non-sparking type, the item itself is not explosion-protected.**
-  **Before connecting your item with the local power supply, make sure the item has been designed for your local voltage; refer to data plate on item.**
-  **Connect your item with a grounded outlet only.**

E

Turn power switch to OFF when item is not in use, or before disconnecting from mains.

Caution ! Hazard of burns when running the hot bath at temperatures exceeding 60°C. Avoid touching hot edge of bowl, evaporator flask and fluid in bowl. Use gloves when replacing evaporator flask.

Never run hot bath dry.

Repair work is limited to skilled personnel so authorized by Heidolph.

Your item needs a solid stand.


Locate the rotary evaporator on a smooth, horizontal table or on the ROTACOOOL reflux cooler. Make sure not to obstruct area around evaporator.

Avoid overpressure in the system. Explosion hazard. For reasons of safety open aerating valve, when working at ambient pressure, in order to avoid excessive overpressure.

SET-UP

1 Item set-up

1a fig. 1

Locate basic item on a smooth and solid surface, e.g. lab table or on ROTACOOOL reflux cooler. Item is set-up horizontally.

Locate hot bath (2) feet in rails of console (1a). Controls facing forward.

All of the four feet of the hot bath feature height adjustment, to prevent tipping. Re-adjust feet as required until hot bath (2) properly is seated on console (1a).

2 Electric wiring

Before connecting with your local power supply, double check to be sure:

your item has been designed for your local voltage as depicted on the data plate (also refer to "Safety Information" section). Refer to "General", if you use a different outlet / plug system;

E

power switches (1b & 2a) for item and hot bath are OFF;

rotary controls (1c & 2b) for rotation and hot bath temperature are set to ZERO (LH stop).

Connect power cord with outlet. Use connectors installed in the rear panel.

Use junction cable to connect hot bath (2) with the item.


fig. 2

X power cord

Y junction cable, hot bath

3 Install condenser support bar, when using G3 and G5eco glassware kits

 view fig. 1,3

G3 and G5eco glassware kits are protected against tipping with condenser support bar.

First secure bar (3) in threaded hole (4) with hex-key screw (5). Tilt bar (3) to fit the condenser (K) to be held.


fig. 3

Fit cross sleeve (6) on bar (3) and secure with clamp screw.

Plug condenser clamp (7) in cross hole of cross sleeve (6) and secure with clamp screw.

4 Install vapor tube

 view fig. 1,4,5

The vapor tube ass'y. (8) is packed with the glassware kit. First remove coupling ring from flange (12), remove cocking spring (13) and PTFE-seal (11), unthread vapor tube coupling ring (9) from drive (4) and remove clamping sleeve (10) by pulling.

Slip coupling ring (9) and clamping sleeve (10) on vapor tube (8), bead pointing forward.


fig. 4

Make sure, small bead on ID of clamping sleeve (10) snaps into groove of vapor tube (8).

Introduce vapor tube ass'y. (8) into drive head from the right and secure on drive head with coupling ring (9). Pull ring hand-tight.

For pulling tight and loosening the vapor tube's coupling ring (8), drive shaft (4) can be locked with button (23).

Slip PTFE seal (11), sealing edge pointing forward, over vapor tube (8), from the left. Location of sealing edge is very important to guarantee proper sealing action under vacuum.


fig. 5

5 Install glassware set G1

 view fig. 1,5,6

Raise lift

Unscrew coupling ring (12) from drive head and remove cocking spring (13). Slide coupling ring (12) over condenser (K) flange. Slide cocking spring (13) over flange as well. Screw condenser (K) on drive head with coupling ring (12), pressing flange against PTFE seal 26 (11).

Cooling water is connected to studs (W) with fitting (H) and threaded cap GL 14 (G).

Collecting flask (N) is secured to condenser (K) with a clamp (O). The clamp (O) is secured by turning the knurled screw clockwise. Attach evaporator flask (P) to vapor tube 26 (8) with clamp (Q).

Inlet pipe (M) is inserted into ground glass fitting at the far end of the condenser (K). Doing this connection, make sure drip ring is located on the PTFE tube as to prevent condensate from flowing back into evaporator flask (P).

Insert supplied inlet tube (S) through cap (E) and washer (D) mounted on inletcock. Secure tube by turning cap clockwise.

Vacuum is connected to fitting (V) with connector (H) and threaded cap GL 14 (G).

Before lowering the lift, locate hot bath centrally underneath evaporation flask (P).


fig. 6 * location of drip ring

P/N	item	name	qty.
514-00000-01	8	vapor tube 26, NS 29/32	1
514-00020-03	8	vapor tube 26 NS 24/40	1
514-00100-00	K	condenser G1	1
514-00110-00	K	condenser G1B, plastic-coated	(1)
514-51000-00	M	inlet pipe	1
514-74000-00	P	evaporator flask 1000 ml, NS 29/32	1
514-74000-05	P	evaporator flask 1000 ml, NS 24/40	1
514-84000-00	N	collecting flask 1000 ml, S 35/20	1
514-84000-02	N	collecting flask 1000 ml, S 35/20 plastic-coated	(1)
515-40000-00	Q	clamp NS 29/32	1
515-40000-01	Q	clamp NS 24/40	1
515-42000-00	O	clamp S 35/20	1
23-09-03-01-24	L	threaded cap GL 18	1
23-09-03-01-27	G	threaded cap GL 14 (for cooling water and vacuum)	3
23-30-01-04-69	E	Threaded cap GL10 red	1
23-08-06-03-26	D	O-ring 3,2x2,5	1
23-30-01-04-72	S	inlet tube (PTFE-) 3,5/4,5x600	1
11-300-005-22	H	hose fitting for threaded cap GL14	3

Use screwdriver to remove cocking spring (13) from condenser flange G1 (K)

6 Install glassware set G3

view fig. 1,5,7

Raise lift.

Unscrew coupling ring (12) from drive head and remove cocking spring (13). Slide coupling ring (12) over condenser flange (K).

Slide cocking spring (13) over the flange as well. Screw condenser (K) on drive head with coupling ring (12), pressing flange against PTFE seal 26 (11).

Arrange condenser clamp (7) on condenser (K) and secure on column (3) with clamp screws. Match column (3) tilt to condenser (K). For more details about installation of condenser holder, refer to page 31.


fig. 7

* location of drip ring

Cooling water is connected with studs (W), hose fitting (H) and threaded cap (G).

Collecting flask (N) is secured on condenser (K) with a clamp (O). The flask clamp (O) is secured in place by turning the knurled-head screw clockwise. Connect evaporator flask (P) with vapor tube 26 (8); use clamp (Q).

The inlet pipe (M) is inserted into the ground glass stud on the left of the condenser. When doing this connection, make sure drip ring is located on the PTFE tube as to prevent condensate from flowing back into evaporator flask (P).

Trim PTFE hose until reaching bottom of the evaporator flask (P).

Insert supplied inlet tube (S) through cap (E) and washer (D) mounted on inletcock. Secure tube by turning cap clockwise.

Install the vacuum fitting (T) into the ground glass stud on top of condenser (K).

Connect vacuum with stud (V), fitting (H) and threaded fitting GL 14 (G).

P/N	item	name	qty.
514-00000-01	8	vapor tube 26, NS 29/32	1
514-00020-03	8	vapor tube 26, NS 24/40	1
514-00300-00	K	condenser G3	1
514-00310-00	K	condenser G3B, plastic-coated	(1)
514-51000-00	M	inlet pipe	1
514-00001-00	T	vacuum fitting, NS 29/32	1
514-74000-00	P	evaporator flask 1000 ml, NS 29/32	1
514-74000-05	P	evaporator flask 1000 ml, NS 24/40	1
514-84000-00	N	collecting flask 1000 ml, S 35/20	1
514-84000-02	N	collecting flask 1000 ml, S 35/20 plastic-coated	(1)
515-40000-00	Q	clamp NS 29/32	1
515-40000-01	Q	clamp NS 24/40	1
515-42000-00	O	clamp S 35/20	1
23-09-03-01-24	L	threaded cap GL 18	2
23-09-03-01-27	G	threaded fitting GL 14 (for cooling water and vacuum)	3
11-300-005-22	H	connector for the threaded fitting GL14	3
23-30-01-04-69	E	threaded cap GL10 red	1
23-08-06-03-26	D	O-ring 3,2x2,5	1
23-30-01-04-72	S	inlet tube (PTFE-) 3,5/4,5x600	1
591-00060-00		condenser holder	1

7 Install glassware set G5eco

view fig. 1,5,8

Raise lift

Unscrew coupling ring (12) from drive head and remove cocking spring (13). Slide coupling ring (12) over flange of cold-trap (K). Slide cocking spring (13) over flange as well.

Screw cold trap shroud (K) on drive head with coupling ring (12), seating flange against PTFE seal 26 (11).

Refer to fig. 1 for installation of condenser bracket.

Plug bracket ass'y. (S) into horizontal cross-clamp bore (6) and secure with clampscrew. Match stand column (3) with glass kit inclination and fit retainer ass'y. (A) over cold trap shroud (K). Secure cold trap shroud (K) with clamp strap being part of retainer ass'y. (A).

Collecting flask (N) is secured to cold-trap shroud with clamp (O). Secure clamp (O) by turning the knurled-head screw clockwise.

Fasten evaporator flask (P) to vapour tube 26 (8) with clamp (Q).

The inlet pipe (M) is inserted into ground glass stud on the left of the cold-trap shroud. Doing this connection, make sure drip ring is located on the PTFE tube as to prevent condensate from flowing back into evaporator flask (P).

Trim PTFE hose until reaching bottom of evaporator flask (P).

Insert supplied inlet tube (S) through cap (E) and washer (D) mounted on inletcock. Secure tube by turning cap clockwise.

Locate gasket (X= silicon with FEP, or X1= silicon) in sealing groove of cold trap shroud (K). Insert centering ring (B) in cold trap shroud (K) and insert cold trap (Y) in cold trap shroud (K).

Connect vacuum source with stud (V), connector (H) and threaded fitting GL 14 (G).


fig. 8

P/N	item	name	qty.
514-00000-01	8	vapor tube 26, NS 29/32	1
514-00020-03	8	vapor tube 26, NS 24/40	1
514-00520-00	K	cold-trap shroud eco	1
514-00520-01	K	cold-trap shroud eco, plastic-coated	(1)
514-00501-00	Y	cold-trap	1
23-30-01-01-39	X	seal G5 (Silicon/FEP)	1
23-30-01-01-88	X1	Seal G5 (Silicon)	1
514-51000-00	M	inlet pipe	1
514-74000-00	P	evaporator flask 1000 ml, NS 29/32	1
514-74000-05	P	evaporator flask 1000 ml, NS 24/40	1
514-84000-00	N	collecting flask 1000 ml, S 35/20	1
514-84000-02	N	collecting flask 1000 ml, S 35/20 plastic-coated	(1)
515-40000-00	Q	clamp NS 29/32	1
515-40000-01	Q	clamp NS 24/40	1
515-42000-00	O	clamp S 35/20	1
23-09-03-01-27	G	threaded fitting GL 14 (for cooling water and vacuum)	1
11-300-005-22	H	connector for threaded fitting GL14	1
591-00061-00	A	retainer ass`y	1
23-30-01-04-88	B	centering ring	1
23-30-01-04-69	E	Threaded cap GL10 red	1
23-08-06-03-26	D	O-ring 3,2x2,5	1
23-30-01-04-72	S	inlet tube (PTFE-) 3,5/4,5x600	1

8 Cooling water fitting on glassware kit (G1, G3)

view fig. 6,7

Cooling water is connected (except G5eco glassware kit) with the threaded glass stud GL 14 of the cooling coil, using barbs (H) and threaded caps GL 14 (G). Use hose with 7 - 8 mm ID and secure with hose clamps.

When selecting parameters of distillation, you need to make reference to cooling water temperature all the time (refer to Chapter "Remarks about distillation setting").

9 Connect vacuum pressure to glassware kit

view fig. 6,7,8

Connect vacuum pressure with threaded glass stud GL 14, using barbs (H) and threaded caps GL 14 (G). Use vacuum hose with an ID of 7 - 8 mm and secure with hose clamp.

While working with vacuum pressure make reference to local safety regulations and setting instructions.

To guarantee proper vacuum tightness, we recommend to coat all ground mating surfaces with industrial quality silicone grease; another option would be use of conventional PTFE inserts.

10 Install Woulff flask / vacuum valve

Woulff flask and vacuum valve are interlinked directly (refer to fig. 9). To this end use the **PTFE-hose** shipped with the Woulff flask (V1 fig. 9). If you are going to install Woulff flask and vacuum valve separately, use the **PTFE adapter** (V2, fig. 10) for connecting the vacuum hose.

Use brackets from your hardware kit to attach the Woulff flask and / or vacuum valve beside the lift column of the basic item with two screws (fig. 10).


fig. 9

fig. 10

11 Vacuum lines

How to route hoses for the individual vacuum systems:


fig. 11

- | | |
|--------------------------------|----------------------|
| 1 condenser on unit | 7 Woulffs' bottle |
| 2 ventilation | 8 vacuum hose |
| 3 pump exhaust | 9 line vacuum valve |
| 4 exhaust condenser (optional) | 10 power supply |
| 5 Vacuum valve | 11 vacuum connection |
| 6 Y hose connector | |


fig. 12

12 Install protective hood, (option P/N 591.00010.00)

This protective hood (25) protects the user from glass fragments, a possible shock wave will escape via the top end. This hood guarantees adequate protection in both, the operational mode and with the lift in upper position.

The protective hood (25) with its structural member is attached to the vertical stop (15) with 2 screws. For opening or closing the protective hood (25) through 90 degrees, use the handle provided. To displace the protective hood (25) loosen clamp screw (16); this way, you can locate the hood over the center of the hot bath; if, for example (using bigger-size evaporator flasks), the hot bath may be displaced to the right and the protective hood (25) arranged in the right place. Secure with clamp screw (16).


fig. 13

13 Install protective shield, (option P/N 591.00020.00)

With the lift in its upper position, a protective shield (17) will protect the operator from fragments, e.g. when replacing a flask still being under vacuum.

The protective shield is hooked-up at the hot bath's metallic frame (2) and clamped with screw (18) (180° turn). Make sure plane of eccentric is pointing up.

Before clamping, position protective shield (17) in a manner that allows to open the protective hood (25) without collision (try out).


fig. 14

CONTROLS AND OPERATION

Before connecting power cord with outlet, make sure that:

Item's voltage and frequency are the same as your local power supply network (refer to data plate on item).

Power switches are OFF ("0"). control knobs are at their lefthand stop, in order to avoid accidental rotation of the evaporator flask.

14 Adjust evaporator flask angle

Evaporator flask angle may be adjusted by turning the rotary drive.

To this end hold condenser (K) of your glassware kit with your right hand, turn grey knob (19) CCW through at least one revolution ("unlock"), carefully tilt condenser (K) to adjust evaporator flask angle.

fig. 15

Pull rotary knob (19) tight by turning clockwise ("lock").

15 Use of hand lift

The hand lift features a big-size handle (20) at the front end of the item. The lift is balanced for the average weight of a standard glassware kit.

fig. 16

Depress lift handle (20) ("push") to unlock lift mechanism. Move lift handle (20) to the left to lower, or to the right to raise the complete drive mechanism with glassware kit and evaporator flask.

Release handle (20) ("push") to lock lift in the present position.

16 Adjust immersion depth of evaporator flask

The lift is used to immerse the evaporator flask in the bath fluid. To use the same immersion depth for all distilling jobs, we recommend to firmly set immersion depth.

To this end loosen black rotary knob (22, fig. 15) on lift column by turning CCW ("unlock"). Vertical stop (15, fig. 17) drops (sustain by pulling down). Descend lift until reaching proper immersion depth in heating bath (2). It may be required to sustain this action by hand. Pull rotary knob tight (22, fig. 15) by turning clockwise ("lock"). At this point, immersion depth of the evaporator flask is set permanently.

fig. 17

If immersion depth is not satisfactory after adjustment, repeat above steps.

 If protective hood is installed, you will need to displace vertical stop (15) by hand:

Loosen vertical stop (15) with black rotary knob (22, fig. 15) ("unlock"), hold protective hood's frame (14, fig. 18) and move vertical stop till hitting upper end. Pull rotary knob (22, fig. 15) tight ("lock").

Descend lift until reaching immersed depth wanted.

Hold protective hood's frame with your hand, the other hand is going to loosen the black rotary knob (22) ("unlock"), and slip vertical lock (15) down to its extreme position. Pull rotary knob (22) tight ("lock").

 Protective hood goes on the vertical stop (15), refer to fig. 17. The lift being lowered will guarantee permanent immersed depth in the hot bath, as specified before.


fig. 18

17 Evaporator flask draw-off fixture

An evaporator flask draw-off system allows to loosen a stuck evaporator flask (P) from the vapor tube (8) at ease.


fig. 19

Remove flask clamp (Q). Hold evaporator flask with your right hand (use gloves to avoid burns). With your left hand turn draw-off nut (24) CCW until evaporator flask is drawn off.

Hold evaporator flask (P), do not drop into hot bath (flask might be broken, hazard of splashes).

Evaporator flask (P) may be hot ! Allow to chill before handling, or use gloves.

Allow hot bath to chill, in order to avoid burns.

18 Displace heating bath

Using bigger-size evaporator flasks or couplings between evaporator flasks and vapor tube (8, fig. 19), you will need to increase clearance between hot bath and drive head.

All of the LABOROTA –models feature feet in **rails**, the heating bath can be displaced as required.


fig. 20

After displacing, re-adjust immersed depth and angle of the evaporator flask.

For extra extension beyond the range of displacement, an optional extension panel (P/N 591-00090-00) is available; it attaches on the item with ease to give an added extension range of displacement to 160 mm.

19 Fill hot bath

Using deionized or distilled water for filling the hot bath, you need to add 0.2 % of borax ($\text{Na}_2\text{B}_4\text{O}_7 \cdot 10\text{H}_2\text{O}$) as corrosion inhibitor.

Please note: The use of tapwater as heating media can lead to calcification which can provoke crevice corrosion. Therefore the user is strongly recommended to clean the heating bath with an adequate polish on a regular basis.

Never run your heating bath (2) dry !

Before heating, fill heating bath (2) with fluid. The bowl features min- and max. marks (2c).

Max. fill level may be reached with the evaporator flask immersed only. Using big-size evaporator flasks carefully immerse item in order to avoid spillage.

For ease of transport, the hot bath (2) features rubber grips (2d) to improve friction.

Hazard of burns when running hot bath (2) at temperatures of more than 60° C. Beware of metal edge of hot bath, fluid and evaporator flask.


fig. 21

20 Work with protective hood and protective shield (optional, P/N 591.00010.00 and 591.00020.00)

If you cannot locate the rotary evaporator behind a glass shield, optional protective hood (25) and shield (17) offer adequate protection for the user.

The protective hood raises together with the lift, hence giving room for replacing the evaporator flask (P).

Folding the protective hood up through 90° with its handle (25a) the evaporator flask's section (P) is accessible. As described in Chapter 12, the protective hood easily adjusts to fit position of the hot bath, e.g. if the hot bath needs to be displaced.

Do not open protective hood in lower position of lift, if you use hood and shield together (collision between hood and shield).

21 Add substance to be distilled

Substance to be distilled may be added through feeder tube (M) under vacuum.

Trim inlet tube (S) (ref. fig. 6-8) and immerse end of hose in container with the substance to be distilled.

Open feed with plug of cock (M) (turn aperture of cock towards nipple). When closing cock (M) make sure not to turn the plug through more than 90 degrees.

22 Aerate by hand

The cock in feeder tube (M) also allows system aeration by hand. To this end turn aperture in cock's plug up towards aerating hole (opposite hose liner).

23 Adjust hot bath temperature

Turn hot bath ON with push-button (2a); button lighted green.

Turn power switch (1b, fig. 1) to ON before starting adjusting work.

Select heating bath temperature with rotary knob (2b).

Temperature marks arranged around rotating knob (2b) just are a rough

reference for setting nominal temperature.

A yellow pilot light (2e) on top of the switch lights when item is heating.


fig. 22

For reasons of safety, hot bath temperature is monitored by a temperature limiter. If temperature control fails and actual temperature exceeds max. temperature by about 10 %, the temperature limiter permanently disconnects heater.

Damage to the guard hood! Set a maximum setpoint of 120 °C for operation with a guard hood.

This temperature limiter is located on the underside. Disconnect status is displayed by a red button (2f) protruding over the base panel by about 2 mm. To reset temperature limiter just depress red button (2f) after some time for chilling off.


fig. 23

The temperature limiter having dropped, is an indication to have the hot bath inspected by a service engineer authorized by Heidolph.

24 Turn on and set rotary speed

First turn rotary evaporator's power switch (1b) to ON.

Before turning power switch to ON, make sure temperature control knob (1c) is in extreme left position to avoid accidental rotation of the evaporator flask.

Adjust rotational speed with knob 1c. Refer to scale for speed setting. In extreme left position of this knob, speed is zero.

INFORMATION ABOUT DISTILLATION SETTINGS

Using a rotary evaporator for thermal cracking of substances some basic consideration should be made to achieve the best results possible. Optimize settings as described below:

25 Rotary speed of evaporator flask

Raise distillation speed by raising rotational speed. This methods cuts distillation time and protects the substances handled.

26 Temperature difference between hot bath and vapor

This temperature difference should at least be 20 K to reach adequate distillation speed. If acceptable by thermal stability of the substances handled, you even may raise temperature difference.

Rule of thumb: double temperature to double distillation speed.

27 Temperature difference between vapor temperature and coolant

This temperature difference should be no less than 20 K to guarantee adequate condensation. Low temperature difference, effectiveness of solvent recovery drops. Choose an as high as possible temperature difference in particular for high-boiling substances.

28 Solvent data

Examples for using reference table and nomogram:

This reference table details substance data needed for distillation, whereas the nomogram reflects relationship between pressure and boil temperature of some solvents. Go by the Clausius-Clapeyron equation and plot $1/T$ against $\log p$.

28.1 At boil temperature of 40°C

At the vacuum controller: in column "vacuum for boil temperature" use the "Set Vac" nominal value as depicted for 40°C.

28.2 At boil temperatures deviating from 40°C

On temperature axis of the nomogram mark boil temperature needed.

Draw a line to the right and determine intersection with the solvent line.

From this intersection draw a straight line down to determine the vacuum pressure needed for this particular application.

28.3 Solvents not depicted

The following items might help in determining the right vacuum pressure:


Rise of the straight line is a function of evaporation enthalpy. It is similar for similar chemical substances with similar boil temperature. Hence, straight lines depicted may be used for rough reference to a lightly deviating boil point.

A water jet pump or diaphragm-type pump may lower boil temperature by about 100°C.

Rule of thumb: decrease pressure to ½ will decrease boil point by about 15°C.

Solvent	Total formula	MW [g/mol]	MW [°C]	ΔH_{vap} [J/g]	Vacuum for Boiling point at 40°C	
					[mbar]	[mm(Hg)]
Aceton	C ₃ H ₆ O	58,08	56,5	550	556	387
Acetonitril	C ₂ H ₃ N	41,05	81,8	833	230	173
Benzol	C ₆ H ₆	78,11	80,1	549	236	177
n-Butanol (Butylalkohol)	C ₄ H ₁₀ O	74,12	117,5	619	25	19
tert.-Butanol (tert.-Butylalkohol)	C ₄ H ₁₀ O	74,12	82,9	588	130	98
2-Butanon (Methylethylketon)	C ₄ H ₈ O	72,11	79,6	473	243	182
tert.-Butylmethylether	C ₅ H ₁₂ O	88,15	55,0			
Chlorbenzol	C ₆ H ₅ Cl	112,60	132,2	375	36	27
Cyclohexan	C ₆ H ₁₂	84,16	80,7	389	235	176
1,2-Dichlorethan	C ₂ H ₄ Cl ₂	98,96	82,4	336	210	158
1,2-Dichlorethylen (cis)	C ₂ H ₂ Cl ₂	96,94	59,0	320	479	134
1,2-Dichlorethylen (trans)	C ₂ H ₂ Cl ₂	96,94	47,8	313	751	563
Dichlormethan (Methylenchlorid)	CH ₂ Cl ₂	84,93	40,7	373	atm.	atm.
Diethylether	C ₄ H ₁₀ O	74,12	34,6	392	atm.	atm.
Diisopropylether	C ₆ H ₁₄ O	102,20	67,5	318	375	281
Dimethylformamid	C ₃ H ₇ NO	73,09	153,0		11	8
1,4-Dioxan	C ₄ H ₈ O ₂	88,11	101,1	406	107	80
Ethanol	C ₂ H ₆ O	46,07	78,4	879	175	131
Ethylacetat	C ₄ H ₈ O ₂	88,11	77,1	394	240	180
Heptan	C ₇ H ₁₆	85,09	98,4	439	120	90
Hexan	C ₆ H ₁₄	86,18	68,7	370	335	251
Methanol	CH ₄ O	32,04	64,7	1225	337	253
3-Methyl-1-Butanol (Isoamylalkohol)	C ₅ H ₁₂ O	88,15	130,6	593	14	11
Pentachlorethan	C ₂ HCl ₅	202,30	160,5	203	13	10
Pentan	C ₅ H ₁₂	72,15	36,1	382	atm.	atm.
n-Pentanol (Amylalkohol)	C ₅ H ₁₂ O	88,15	137,8	593	11	8
1-Propanol (n-Propylalkohol)	C ₃ H ₈ O	60,10	97,8	787	67	50
2-Propanol (Isopropylalkohol)	C ₃ H ₈ O	60,10	82,5	701	137	103
1,1,2,2-Tetrachlorethan	C ₂ H ₂ Cl ₄	167,90	145,9	247	35	26
Tetrachlorethylen	C ₂ Cl ₄	165,80	120,8	233	53	40
Tetrachlormethan (Carbontetrachlorid)	CCl ₄	153,80	76,7	225	271	203
Tetrahydrofuran	C ₄ H ₈ O	72,11	66,0		357	268
Toluol	C ₇ H ₈	92,14	110,6	425	77	58
1,1,1-Trichlorethan	C ₂ H ₃ Cl ₃	133,40	74,1	251	300	225
Trichlorethylen	C ₂ HCl ₃	131,40	86,7	265	183	137
Trichlormethan (Chloroform)	CHCl ₃	119,40	61,3	263	474	356
Wasser	H ₂ O	18,02	100,0	2259	72	54
Xylol (Isomeren-Gemisch)	C ₈ H ₁₀	106,20	137-143	390	25	19

Conversion factor [mbar] to [mm(Hg)]: $[mm(Hg)] \approx \frac{3}{4}[mbar]$

CLEANING AND MAINTENANCE

For cleaning, wipe housing with a damp cloth (add some sort of mild liquid soap).

NOTE: To avoid damage to the surface finish, avoid using chlorine bleach, chlorine-based cleaners, abrasive substances, ammonia, rags or cleaning pads containing metallic particles.

Please note: The use of tapwater as heating media can lead to calcification which can provoke crevice corrosion. Therefore the user is strongly recommended to clean the heating bath with an adequate polish on a regular basis.

Vapor tube (8, fig. 5) and PTFE seal (11, fig. 5) need regular upkeep. To this end, remove, clean, and inspect vapor tube and PTFE seal (sealing edge); replace PTFE seal on condition.

The item is maintenance-free. Repair work is limited to engineers authorized to do so by HEIDOLPH. So therefore get in touch with your local HEIDOLPH dealer or field representative.

DISASSEMBLE, FORWARD & STORE

Disassemble

Turn the item OFF and disconnect mains plug.

Transport and store

Use original box or a similar pack to store your item in order to prevent damage in transit. Tape box as required.

Store item in a dry place.

CAUTION: avoid shocks and vibration while in transit.

DISPOSAL

For disposal, please comply with your local or national regulations.

Split by metal, plastic, etc.

Packing material to be treated as described above (material split).

TROUBLESHOOTING

No heat:

- check power cord
- power switch (1b) broken (call service engineer)
- max. temperature limiter dropped (refer to Chapter "Adjust hot bath temperature")
- hot bath heater broken (call service engineer)
- LABOROTA master switch (1b) OFF; (turn master switch 1b ON)

Drive won't run:

- check power cord
- master switch (1b) broken (call service engineer)
- speed control knob at lefthand end, turn right
- motor problem (call service engineer)

Insufficient vacuum:

- check seals and fittings
- check ground fittings, lubricate on condition
- check PTFE seal for proper installation (sealing edge faces drive head). Refer to "set-up, vapor tube and seal",
- check vacuum pump's suction line

SPECIFICATIONS

model	LABOROTA 4000eco
drive	condenser motor with electronic speed control
speed range	20-270/min
speed control	scale
lift actuation	manual
lift speed	manual
vertical travel	140 mm
heating power	1300 W
temperature range, hot bath	20-100°C
accuracy, hot bath	±2.5 K
heating bath temperature control	analog
material, hot bath	V4A (1.4404)
heating bath diameter	240 mm
power supply	230V/50/60Hz oder 115V/50/60Hz
dissipated power	1400 W
weight, glassware not included	17 kg
size (wxdxh), lift in position low, glassware not installed	396 x 440 x 360 mm space requirement, includes lift lever
power cord	2 m
protective class	IP 20
ambient conditions	0-40°C at 80% rel. humidity
evaporative power	abt. 1000 ml H ₂ O/h
cooling area	1200 cm ²

WARRANTY, LIABILITY & COPYRIGHT

Warranty

Heidolph Instruments warrants that the present product shall be free from defects in material (except wear parts) and workmanship for 3 years from the date shipped from the manufacturer's warehouse.

Damage in transit is excluded from this warranty.

To file for such warranty service, contact Heidolph Instruments (phone ++49-9122-9920-68) or your local Heidolph Instruments Dealer. If defects in material or workmanship are found, your item will be repaired or replaced at no charge.

Misuse, abuse, neglect or improper installation are not covered by this warranty. Alterations to the present warranty need Heidolph Instruments' consent in writing.

Exclusion Clause

Heidolph Instruments cannot be held liable for damage from improper use or misuse. Remedy for consequential damage is excluded.

Copyright

Copyright in pictures and wording of the present Instruction Manual is held by Heidolph Instruments.

FAQ / REPAIR WORK

If any **aspect** of installation, operation or maintenance remains unanswered in the present Manual, please contact the following address:

For **repair** services please call Heidolph Instruments (phone: +49 - 9122 - 9920-68) or your local, authorized Heidolph Instruments Dealer.

**Note**

You will receive approval for sending your defective item to the following address:

Heidolph Instruments GmbH & Co. KG

Lab Equipment Sales

Walpersdorfer Str. 12

D-91126 Schwabach / Germany

phone: ++49-9122-9920-68

Fax: ++49-9122-9920-65

E-Mail: Sales@Heidolph.de

**Safety Information**

When shipping items for repair that may have been contaminated by hazardous substances, please:

Advise exact *substance*.

Take proper *protective* measures to ensure the safety of our receiving and service personnel.

Mark the pack IAW Hazardous Materials Act.

CE-DECLARATION OF CONFORMITY

We herewith declare that the present product complies with the following standards and harmonized documents:

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© Heidolph Instruments GmbH & Co. KG