Temperature

Controlled Shaking

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Operating instructions

Hei-MIX Unimax 1010
Promax 1020
Duomax 1030
Polymax 1040

 **heidolph**
research made easy

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Seite 4 – 27

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Translation of the original instructions
Page 32 – 55

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Contents

Introduction

About this document	32
Typographic conventions	32
Copyright protection.....	32

General notes

Basic product information	33
Guidelines applied, product certification	33
California Residents	33
Residual risk.....	33
Intended use	33
Reasonably foreseeable misuse.....	33
Transportation.....	34
Storage.....	34
Acclimatization.....	34
Permissible ambient conditions.....	34

Safety

General safety information.....	35
Electrical safety.....	35
Operational safety	35
Work safety	36
Personal protective equipment (PPE).....	36
Environmental protection.....	36
Biohazard.....	36
Special hygiene measures for the use of laboratory equipment in food, cosmetics and pharmaceutical production	37
General Measures	37
Device-specific measures	37
Other regulations	38

Device description

Mechanical design	39
-------------------------	----

Commissioning

Set up the device	40
Power supply	40
Attachments	41

Operation

Switch the device on/off	42
Unimax 1010, Promax 1020	42
Duomax 1030, Polymax 1040	42
Operating modes	43
Continuous operation Unimax 1010, Promax 1020.....	44
Timer mode Unimax 1010, Promax 1020.....	45
Continuous operation Duomax 1030, Polymax 1040.....	46
Timer mode Duomax 1030, Polymax 1040.....	46

Troubleshooting

Troubleshooting.....	47
----------------------	----

Annex

Technical Specifications	48
Performance range/admissible load	50
Scope of delivery	51
Accessories.....	51
Service	52
General cleaning instructions	52
Replace the device fuse	52
Repairs	53
Maintenance.....	53
Disposal.....	53
Contact information Heidolph international.....	54
Warranty Statement	54
Certificate of decontamination	55

About this document

This operating manual describes the features and operation of platform shakers type Hei-MIX Unimax 1010, Hei-MIX Promax 1020, Hei-MIX Duomax 1030 and Hei-MIX Polymax 1040. The operating instructions manual is an integral part of the delivery!

Typographic conventions

Standardized symbols, highlighting elements, and signal words are used in this document to identify warnings, cautions, important information, and special text contents.

Symbol	Signal word / explanatory note
	<p>Warning symbols in combination with a signal word indicate dangers:</p> <p>DANGER Indicates a hazardous situation which, if not avoided, will result in death or serious injury.</p> <p>WARNING Indicates a possible hazardous situation which, if not avoided, may result in death or serious injury.</p> <p>CAUTION Indicates a hazardous situation which, if not avoided, could result in minor or moderate injury or material/environmental damage.</p>
	<p>Mandatory signs are used to indicate important information regarding the product handling.</p> <p>This information is used to ensure operational safety and to maintain the value of the product.</p>
[GUI]	<p>Parameter designations, display texts, and device labels are highlighted in text and tables in a typographic manner to facilitate the assignment on the device.</p>
→	<p>The arrow symbol indicates instructions to be followed in order to ensure the operational safety when handling the product.</p>

Copyright protection

This publication is protected by copyright and intended for internal use by the purchaser of the product only.

No part of this publication may be transmitted or reproduced in any form, by any means, without the prior written consent of the copyright owner Heidolph Instruments GmbH & Co. KG. Any violation is subject to compensation for damages.

Basic product information

Guidelines applied, product certification

CE Marking

The device complies with the following standards:

- Machinery Directive 2006/42/EC
- EMC Directive 2014/30/EU

California Residents

Important information for California residents regarding Prop 65. Please visit www.P65Warnings.ca.gov for more information.

Residual risk

The device was designed and manufactured in accordance with the latest technical standards at the time of development and the recognized safety regulations. During installation and use, as well as during maintenance work, repairs and cleaning, there are nevertheless certain residual risks associated with the device described.

These are identified and described at the appropriate points in this document.

Intended use

The platform shakers described in this document have been specially developed for the following tasks:

- Shaking
- Mixing
- Emulsifying
- Suspending
- Separating
- Dissolving
- Staining

The areas of application of the platform shakers described include chemical, biological and environmental-analytical laboratory and research applications, applications of basic research and comparable facilities.

Due to its design, the device in its delivery condition may only be used in analytical processes or in laboratory-like conditions in the food, cosmetics, and pharmaceutical industries as well as other comparable industries that manufacture products intended for consumption by humans or animals, or for use on humans or animals.

Any other use of this device is not considered as intended!

Reasonably foreseeable misuse

For use under conditions or for purposes deviating from the intended use, additional measures may become necessary, and/or specific guidelines and safety regulations will have to be observed (see section "Special hygiene measures for the use of laboratory equipment in food, cosmetics and pharmaceutical production" on page 37). Corresponding requirements must be evaluated and observed by the operator in each individual case.

Compliance with and implementation of all relevant guidelines and safety measures for the respective field of application is within the sole responsibility of the operator.

All risks resulting from improper use are solely borne by the operator.

The device may exclusively be operated by authorized and instructed personnel. Training and qualification of the operating personnel as well as ensuring that the device is operated with responsibility are the sole responsibility of the operator!

Transportation

During transport, avoid severe shocks and mechanical stresses that can cause damage to the device.

Keep the original packaging in a dry and protected place for later use.

Storage

Always store the device in its original packaging. To protect against damage and unreasonable material aging, store the device in a dry environment that should be as temperature-stable and dust-free as possible.

Recommended ambient conditions for storage:

- 5 °C – 31 °C, up to 80 % rel. humidity
- 32 °C–40 °C, decreasing linearly up to 50 % rel. humidity

Acclimatization

After each transport and after storage under critical climatic conditions (e.g. high temperature difference between inside and outside), allow the device to acclimatize at room temperature for a minimum of two hours to prevent possible damage from condensation before putting it into operation at the place of use. If necessary, extend the acclimatization phase if the temperature differences are very high.

Make all supply connections (power supply, tubing) only after the device has been acclimatized!

Permissible ambient conditions

The device is designed for indoor use only. Permissible ambient conditions for operation:

- 5 °C – 31 °C, up to 80 % rel. humidity
- 32 °C–40 °C, decreasing linearly up to 50 % rel. humidity
- Maximum height above sea level: 2,000 m

When used in corrosive atmospheres, the service life of the device may be reduced depending on the concentration, duration and frequency of exposure.

The device **IS NOT** suitable for outdoor use!

The device **IS NOT** suitable for use in hazardous areas!

General safety information

- Before commissioning and using the device, familiarize yourself with all the safety regulations and guidelines for occupational safety applicable at the place of use and observe them at all times.
- Only operate the device if it is in perfect technical condition. In particular, make sure that the device is free of visible damage.
- If there is missing or misleading information on the device or on occupational safety, contact the responsible safety specialist or our technical service.
- Only use the device in accordance with the regulations for intended use ("Intended use" on page 33).

Electrical safety

- Ensure that the voltage indicated on the rating plate matches the supply voltage of the country in which the device is being used.
- Ensure that the mains socket-outlet is protected by means of a residual-current device (RCD).
- Always use the supplied power supply 3-pin cord provided with the device (phase, neutral, ground).
- Only replace the device fuses with original fuses from the manufacturer (available as accessories, see section "Accessories" on page 51).
- Before each operation, make sure that the power supply cord is free of visible damage.
- Have repairs and/or maintenance work on the device carried out exclusively by an authorized and skilled electrician or by the technical service department of Heidolph Instruments.
- Always switch the device OFF and disconnect it from the power supply before carrying out maintenance work, cleaning, or repairs.

Operational safety

- Operate the device under a closed ventilated fume hood when working with potentially hazardous substances (according to EN 14175 and DIN 12924).
- Do not make any unauthorized changes or modifications to the device!
- Only use genuine spare parts and accessories, or those expressly approved by the manufacturer!
- Rectify malfunctions or faults on the device immediately.
- Shut down the device and disconnect it from the power supply if it is not possible to eliminate the malfunction or rectify the fault immediately.
- Observe all other applicable regulations such as laboratory and workplace guidelines, recognized safety technology rules and special local regulations.

Work safety

- Always use the prescribed personal protective equipment (PPE) such as protective clothing, safety goggles, protective gloves, safety shoes, etc.
- Do not operate any other devices in the immediate vicinity of the device ...
 - which can generate electromagnetic fields in the frequency range between 9×10^3 Hz to 3×10^{11} Hz,
 - which generate emission or radiation sources in the frequency range 3×10^{11} Hz to 3×10^{15} Hz (in the optical spectral range wavelengths from 1,000 μm to 0,1 μm),
 - which generate ultrasonic or ionizing waves.
- Do not operate the unit when adiabatic compression or shock waves may occur (pressure wave ignition).
- Do not process substances that could release energy in an uncontrolled manner (exothermic reaction, spontaneous ignition).

Rectify malfunctions or faults on the device immediately.
 Shut down the device and disconnect it from the power supply if it is not possible to eliminate the disturbance or rectify the fault.
 In such a case, contact the responsible technical service department.

Personal protective equipment (PPE)

The operator must determine and provide the necessary PPE, depending on the respective application and the media and chemicals used.

The corresponding instruction of the personnel is solely within the operator's responsibility.

Environmental protection

When processing environmentally hazardous substances, take appropriate measures to avoid risks to the environment.

The evaluation of corresponding measures such as the marking of a hazardous area, their implementation, and the training of the responsible personnel is the sole responsibility of the operator!

Biohazard

When processing biohazardous substances, take appropriate measures to prevent hazards to persons and the environment, including:

- Instruction of the personnel regarding the necessary safety measures.
- Provision of personal protective equipment (PPE) and instruction of the personnel in its use.
- Marking of the device with a biohazard warning symbol.

The evaluation of corresponding measures such as the marking of a hazardous area, their implementation, and the training of the responsible personnel is the sole responsibility of the operator!

Special hygiene measures for the use of laboratory equipment in food, cosmetics and pharmaceutical production

When laboratory equipment is used in the production processes of the food, cosmetics or pharmaceutical industry, special hygiene measures must be taken by the user to avoid sample contamination and to minimize any risk to humans and the environment as far as possible.

Please observe the following recommendations:

General Measures

- Ensure a clean working and storage environment when handling substances and materials.
- Train all employees in the field of occupational hygiene, document all training measures and check the implementation of all required hygiene measures during operation regularly.
- Use a hygiene control concept such as HACCP (Hazard Analysis and critical Control points). The HACCP comprises the following criteria:
 - Hazard analysis
 - Identification of critical control points
 - Definition of critical limit values
 - Establishment of a system for monitoring and controlling critical hazard control points (CCP)
 - Corrective actions for uncontrollable CCP
 - Establishment of a system to verify the implementation of all HACCP measures
 - Establishment of a system for documenting all associated procedures and protocols

The evaluation of the applicability of the mentioned rules and regulations is within the sole responsibility of the operator!

Device-specific measures

- Regularly clean components that come into contact with the product, such as flasks, seals, tubes, etc. in the autoclave (if available or possible) or chemically (e.g. with ethanol) to sterilize all surfaces.
- Make sure that even products that are intended for single use only are of sufficient purity.
- Do not use open containers.
- Avoid contamination by handling contaminated vessels, apparatus or aids with care.

Contact information

For further information, please contact our after sales service at any time.

Phone: +49-9122-9920-0

Mail: sales@heidolph.de

Other regulations

In addition to the notes and instructions in this document, observe all other applicable regulations such as laboratory and workplace guidelines, hazardous substances ordinances, recognized rules of safety engineering and occupational medicine as well as particular local regulations!

Noncompliance will invalidate any warranty against Heidolph Instruments.


The operator is solely liable for all damage resulting from unauthorized changes or modifications to the unit, from the use of unauthorized or non-genuine spare parts and accessories, or from disregarding the safety instructions and hazard warnings or the manufacturer's instructions!

Mechanical design

From a design point of view, the Hei-MIX Unimax 1010 is identical to the Hei-MIX Promax 1020. However, both devices offer different motion patterns (Hei-MIX Unimax 1010: orbital movement; Hei-MIX Promax 1020: reciprocating movement).

From a design point of view, the Hei-MIX Duomax 1030 is identical to the Hei-MIX Polymax 1040. However, both devices offer different motion patterns (Hei-MIX Duomax 1030: rocking movement; Hei-MIX Polymax 1040: waving movement).

Hei-MIX Unimax 1010
Hei-MIX Promax 1020

Hei-MIX Duomax 1010
Hei-MIX Polymax 1020

Set up the device

CAUTION

Material damage, loss of production

Improper installation of the device may result in direct and indirect damage to property due to falling and/or spilling of fluids!

- Observe the instructions for the correct positioning of the device.
- Observe the generally valid rules for the safe installation of work equipment.

- Place the device in a suitable workplace.
- The working surface must be even and flat.
- The device must be freely accessible at all times.
- Ensure that a sufficient safety distance can be maintained to moving and/or hot components.

The professional mounting and positioning of the device including the supplied accessories is within the sole responsibility of the operator!

Heidolph instruments accepts no responsibility for direct and/or indirect personal injury or damage to property resulting from failure to observe the instructions for the correct mounting and positioning of the device.

Power supply

DANGER

Electric Shock

The device may only be supplied with the required mains voltage via a grounded mains socket-outlet (voltage see rating plate).

IEC appliance inlet with fuse holder

The IEC appliance inlet and fuse holder are located on the back of the device (example: Unimax 1010):

Connect the power supply cord

- Before connecting the power supply cord, make sure that the main switch of the device is in position **0** (off).
- Connect the cable coupling of the supplied coded power supply cord (see figure) to the IEC appliance inlet on the back of the device:

- Connect the power supply cord to a properly secured mains socket-outlet.

Disconnect the power supply cord

- Before disconnecting the mains connection cord, switch off the device.
- Disconnect the mains connection cord from the mains socket-outlet.
- Disconnect the power supply cord from the IEC appliance inlet on the back of the device.

Attachments

Use the various optional attachments to fix your laboratory vessels. When mounting the attachment on the shaking platform of the device, refer to the mounting instructions delivered with each attachment.

At low shaking frequencies, vessels such as Petri dishes can also be placed directly on the non-slip contoured rubber mat.

- Always place individual vessels in the center of the shaking platform.
- Always distribute several vessels evenly on the shaking platform.

Switch the device on/off


Unimax 1010, Promax 1020

- Set the speed control to zero before switching on the device!
- Use the main switch to switch the device on or off, see also section "Mechanical design" on page 39.
 - When on, the main switch lights up green.

Duomax 1030, Polymax 1040

- Set the speed control to zero before switching on the device!
- Use the main switch (located on the right side of the control panel) to turn the device on or off, see also section "Mechanical design" on page 39.
 - When on, the main switch lights up green.

Operating modes

Regardless of the operating mode, observe the following specific safety instructions.

WARNING

Risk of injury, risk of damage to property due to uneven load, overload, leaking fluids and/or falling vessels

- When possible, use closed vessels for processing corrosive, toxic or biohazard substances and seal them safely.
- Increase the speed only gradually to the desired shaking frequency, especially with open vessels, and observe the fluid movements.
- Observe the maximum load capacity of the device platform, see section "Performance range/admissible load" on page 50.
- Lower the shaking frequency or reduce the total load on the shaking platform when the device starts to move during operation.
- Always ensure that the vessels are evenly distributed on the shaking platform.
- Pay particular attention to an even distribution of the weight load when vessels of different sizes and/or differently filled vessels are placed on the shaking platform at the same time!
- Before switching on the device, make sure that all vessels are properly fixed on the shaking platform.
- Use suitable attachments for the vessels used (see also section "Accessories" on page 51).
- Always use the necessary and appropriate personal protective equipment!

Automatic restart

When the mains is restored after a power interruption, the device automatically starts up again with the set parameters!

Depending on the setting and structure, there is a risk that the device will move uncontrolled and/or fluid will be sprayed out of the vessels.

- In the event of an interruption to the mains supply, switch off the device via the main switch.

All devices can be operated in continuous operation or time-controlled via the integrated timer (timer operation).


- The Unimax 1010 and Promax 1020 devices work with a digital timer.
- The Duomax 1030 and Polymax 1040 devices work with a mechanical timer.

The individual functions are described in detail in the following sections.

Continuous operation Unimax 1010, Promax 1020

In continuous operation, the shaking movement can be started and stopped manually as required after switching on the device.

- Place the vessel(s) on the platform. Observe the instructions in section "Attachments" on page 41.
- Ensure that the speed control is set to zero.
- Switch on the device (see "Switch the device on/off" on page 42).
- Press the two timer buttons (up/down) at the same time to set the timer to **000**:

- Press the start/stop button to activate the standby mode.
 - The **start/stop** and **Speed** LEDs light up.
- Start shaking operation by increasing the intensity up to the desired value, using the speed control.
 - The shaking intensity (speed) is shown in the display.
- Stop the shaking movement by turning the speed control back to zero.

The shaking intensity (speed) can be adjusted at any time during operation.

Timer mode Unimax 1010, Promax 1020

In timer operation, a specific period of time can be set for continuous shaking.

- Place the vessel(s) on the platform. Observe the instructions in section "Attachments" on page 41.
- Ensure that the speed controller is set to zero.
- Switch on the device (see "Switch the device on/off" on page 42).
- Use the two timer buttons to set a value in the range of 0 – 999 minutes (LED **Time** flashes).

Press and release the buttons to increase or decrease the display value in one-step increments.

Press and hold the buttons to increase or decrease the display value in fast forward.

The set timer value is retained until the device is switched off and can be adjusted at any time during operation.

- Press the start/stop button to activate the motion mode.
 - The **start/stop** and **Time** LEDs light up.
 - The the time already elapsed is displayed on the device.
- Adjust the shaking intensity with the speed control as required.
 - The shaking intensity (speed) is shown in the display.

In case of stopping the movement by pressing the start/stop button in timer mode, the timer is stopped but not reset. The timer continues to run accordingly when the shaking movement is restarted, starting from the last count value.

In the event of a power interruption or when the device is switched off via the main switch, the timer is reset to zero!

- Use the Speed/Time mode selector to toggle between the timer and the speed display (LED Time/Speed) while the system is running.
- Press Stop/stop button again to deactivate the motion mode.

Continuous operation Duomax 1030, Polymax 1040

In continuous operation, the shaking movement can be started and stopped manually as required after switching on the device.

- Place the vessel(s) on the platform. Observe the instructions in section "Attachments" on page 41.
- Ensure that the speed control is set to zero.
- Switch on the device (see "Switch the device on/off" on page 42).
- Start shaking operation by increasing the intensity up to the desired value, using the speed control.
- Stop the shaking movement by turning the speed control back to zero.


The shaking intensity (speed) can be adjusted at any time during operation.

Timer mode Duomax 1030, Polymax 1040

In timer operation, a specific period of time can be set for continuous shaking.

- Place the vessel(s) on the platform. Observe the instructions in section "Attachments" on page 41.
- Switch the device to time mode, using the time/speed selection switch:

- Use the **Time control** to set the desired operating time.
- Switch on the device (see "Switch the device on/off" on page 42)..
- Start shaking operation by increasing the intensity up to the desired value, using the speed control.

The shaking movement stops automatically as soon as the timer has expired.
The shaking movement can be stopped manually at any time by resetting the speed control to zero.

Troubleshooting

Failure	Cause	Possible remedy
LED of the on/off switch remains off when switching on	No mains voltage	Check the power supply cord for correct connection Check RCD Check fuses of the device
	LED defective	Contact Heidolph Sales or Service
	Timer expired	Restart the timer if necessary
Shaking movement stops	Motor overheat protection has tripped	Wait approx. 20 minutes and reduce the load on the shaking platform
	Electrical fault (no engine noise)	Contact Heidolph Sales or Service
	Mechanical defect (engine noise audible)	Contact Heidolph Sales or Service

In case of recurring errors please contact the responsible sales department or our technical service. Contact address see "Contact information Heidolph international" on page 54.

Technical Specifications

Model-Specific data

Hei-MIX Unimax 1010

Dimensions (W× H× D)	320 × 125 × 375 mm
Weight	approx. 8 kg
Permissible load	max. 5 kg
Usable area (W × D)	290 × 258 mm
Motion	orbital
Orbit	10 mm
Speed range	30 – 500 rpm

Hei-MIX Promax 1020

Dimensions (W× H× D)	320 × 125 × 375 mm
Weight	approx. 8 kg
Permissible load	max. 5 kg
Usable area (W × D)	290 × 258 mm
Motion	reciprocating
Orbit	32 mm
Speed range	30 – 250 rpm

Hei-MIX Duomax 1030

Dimensions (W× H× D)	320 × 185 × 375 mm
Weight	approx. 8 kg
Permissible load	max. 5 kg
Usable area (W × D)	290 × 258 mm
Motion	rocking, 5/10°
Speed range	2 – 50 rpm

Hei-MIX Polymax 1040

Dimensions (W× H× D)	320 × 195 × 375 mm
Weight	approx. 8 kg
Permissible load	max. 5 kg
Usable area (W × D)	290 × 258 mm
Motion	waving, 5/10°
Speed range	2 – 50 rpm

General device data

Drive	Unimax 1010	EC motor
	Promax 1020	
	Duomax 1030	Asynchronous motor
	Polymax 1040	
Protection class	IP30	
Acoustic pressure	< 50 dB (A), in accordance with IEC 61010	
Engine overheat protection	Unimax 1010	n.a.
	Promax 1020	
	Duomax 1030	self-resetting
	Polymax 1040	

Electrical data


Rated voltage	Unimax 1010	100 – 240 V, 50/60 Hz
	Promax 1020	
	Duomax 1030	230 V, 50 Hz or 115 V, 60 Hz
	Polymax 1040	
Fuse	M 1,25 A (230 V) or T 2,0 A (115 V)	
Protection class	I Ⓢ	
Power input	Unimax 1010	50 W
	Promax 1020	
	Duomax 1030	115 W
	Polymax 1040	

Permissible ambient conditions


Operating and storage temperature	5 °C – 31 °C up to 80 % rel. humidity 32 °C – 40 °C up to 50 % rel. humidity (decreasing linearly)	
Maximum height above sea level	2,000 m	

Performance range/admissible load


Hei-MIX Unimax 1010

Hei-MIX Promax 1020

Hei-MIX Duomax 1030 – Polymax 1040

Scope of delivery

Item	Quantity	Product no..
Hei-MIX Unimax 1010*	1	543-12310-00
Hei-MIX Promax 1020*	1	543-22332-00
Hei-MIX Duomax 1030, 5° inclination*	1	543-32205-00
Hei-MIX Duomax 1030, 10° inclination*	1	543-32210-00
Hei-MIX Polymax 1040, 5° inclination*	1	543-42205-00
Hei-MIX Polymax 1040, 10° inclination*	1	543-42210-00
Power supply cord	1	country-specific
Operating instructions	1	01-005-005-67
Guarantee registration / Certificate of decontamination	1	01-006-002-78
EC Declaration of Conformity	1	01-001-025-06

* the product numbers shown are for 230 V devices for the European market. For product numbers of variants, please contact Heidolph instruments.

Accessories

Item	Quantity	Product no.
Device fuse for 230 V version	1	14-002-015-23
Device fuse for 115 V version	1	14-002-015-45

Further information on available accessories can be found on our website at www.heidolph.com!

Service

When carrying out service work on the device (cleaning, maintenance, repair), observe the general instructions and safety information described in this section.

DANGER

Electric Shock

Live components are installed inside the device.

When opening the device, there is a risk of touching live components.

→ Switch the device's main switch off and disconnect it from the power supply before carrying out maintenance work, cleaning, or repairs.

Penetrating liquid poses the danger of an electric shock.

→ When cleaning, avoid the penetration of liquids.

General cleaning instructions

Wipe all surfaces and the control panel with a damp cloth if necessary. Persistent contamination can be removed with mild soapy water.

CAUTION

Damage to the device

Improper cleaning can damage the surfaces of the device.

Penetrating liquid can damage the electronic components inside the device.

→ Clean the device's surfaces with a soft, lint-free and only slightly moistened cloth.

→ Never use any aggressive or abrasive cleaning agents or aids.

Replace the device fuse

The fuse holder for the two device fuses is located on the rear of the device below the IEC appliance inlet, see section "Mechanical design" on page 39.

DANGER

Electric Shock

Before replacing the fuses, switch off the device and disconnect the power supply cord.

Always replace the two device fuses in pairs with original manufacturer fuses, see section "Accessories" on page 51.

After fuse replacement, check the device for a safe condition according to IEC 61010-1.

Repairs

Repairs to the device may only be carried out by authorized skilled experts!

Unauthorized repairs during the warranty period will result in the loss of the warranty claim.

The owner is solely liable for damage caused by unauthorized repairs.

In case of repair contact an authorized dealer or our technical service, see "Contact information Heidolph international" on page 54.

Include the completed Certificate of decontamination with every device return, see "Certificate of decontamination" on page 55.

Maintenance

The device contains no user-serviceable components. If necessary, in the event of abnormal operating behavior such as excessive noise or heat generation, for example, contact our technical service, see "Contact information Heidolph international" on page 54.

Disposal

-
- When disposing of the device, observe the provisions of the WEEE Directive 2012/19/EU and its transposition into national law in the country of use.
 - When disposing of portable batteries, observe the provisions of the European Battery Directive 2013/56/EU and their transposition into national law in the country of use.
 - Check the device and all components for residues of substances that are hazardous to health, the environment, and biohazardous before disposing.
 - Properly remove and dispose residues of substances that are hazardous to health, the environment and biohazardous!
-

Contact information Heidolph international

Heidolph Instruments North America

Phone: 1-866-650-9604
E-mail: service@heidolph.com
www.heidolphNA.com

Heidolph Instruments United Kingdom

Phone: 01799 - 5133-20
E-mail: service@radleys.co.uk
www.heidolph-instruments.co.uk

Local distributors

To find your local distributor please visit www.heidolph.com

Warranty Statement

Heidolph Instruments assumes a three-year warranty against material and manufacturing defects.

Excluded from the warranty are glass and wear parts, transport damage, and damage resulting from improper handling or non-intended use of the product.

For registered products, the warranty period begins on the date of purchase. Register the product with the enclosed warranty card or on our homepage www.heidolph.com.

For non-registered products, the warranty period begins with the date of the serial production (to be determined by the serial number).

In the event of material or manufacturing defects, the product will either be repaired or replaced free of charge within the warranty period.

Certificate of decontamination

Enclose the Certificate of decontamination, duly completed, with your device return. Submissions without a Certificate of decontamination cannot be processed!

Certificate of decontamination

IN CASE OF RETURNS

Please fill in the required fields.

Note: The sender must package the goods properly and appropriately for transport.

Heidolph Instruments GmbH & Co. KG
Walpersdorfer Straße 12
91126 Schwabach, Germany
Phone: +49 (0) 9122 9920-380
Fax: +49 (0) 9122 9920-19
Email: service@heidolph.de

SENDER

Name _____ First name _____
Company/Institution _____ Department _____
_____ Workgroup _____
Address _____
ZC/City _____
Country _____ Phone _____
Email _____

DEVICE DETAILS

Article number _____ Serial no. _____
Ticket number _____
Reason for sending in _____

Has the device been cleaned, decontaminated/disinfected? Yes No (Please mark as applicable)

If yes, which measures were carried out?

Does this device pose a risk to people and/or the environment due to the processing of substances that are hazardous to health, the environment and/or are biohazardous? Yes No (Please mark as applicable)

If yes, with which substances did the device come into contact?

LEGALLY BINDING DECLARATION

The principal/consignor is aware that they are liable to the agent/consignee for losses or damage incurred due to incomplete and incorrect information.

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Zertifikate – Certifications

EU Declaration of Conformity	84
UKCA Declaration of Conformity	85
RoHS Declaration of Conformity.....	86
China RoHS Certification.....	87

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hereby declare, that the product (from serial number on 200066855) designated below is in compliance with the basic requirements of all applicable EU-directives stated below with regard to design, type of model sold and manufactured by us. This certificate will be invalid if the product is modified without the prior written consent and agreement of the manufacturer.

Reax 2	541-21001-XX		
Reax top	541-10000-XX		
Reax control	541-11000-XX		
Reax 20/4	541-20004-XX		
Reax 20/8	541-20008-XX		
Reax 20/12	541-20012-XX		
Titramax 100	544-11100-XX, 544-11200-XX		
Duomax 1030	543-32105-XX, 543-32205-XX, 543-32210-XX		
Polymax 1040	543-42105-XX, 543-42110-XX, 543-42205-XX, 543-42210-XX		
Polymax 2040	542-40005-XX, 542-40010-XX		
Multi Reax	545-10000-XX		
Vibramax 100	544-21200-XX	Titramax 101	544-11300-XX
Rotamax 120	544-41200-XX	Vibramax 110	544-31200-XX
Titramax 1000	544-12200-XX	Inkubator 1000	549-90010-XX
Unimax 1010	543-12310-XX	Promax 1020	543-22332-XX
Unimax 2010	542-10020-XX	Promax 2020	542-20020-XX

Maschinenrichtlinie / Machinery Directive 2006/42/EG
EMV-Richtlinie / Electromagnetic Compatibility Directive 2014/30/EU
Delegierte (RoHS-) Richtlinie / Delegated (RoHS) Directive 2015/863/EU
Angewandte (harmonisierte) Normen / (Harmonized) Standards applied:
EN ISO 12100:2010, EN 61326-1:2013, EN 61010-1:2010, EN IEC 63000:2018

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UK CA

Declaration of Conformity

In accordance with UK Government guidance

This declaration is issued under the sole responsibility of the manufacturer,
Heidolph Instruments GmbH & Co. KG
Walpersdorfer Straße 12
91126 Schwabach / Germany

Product: Shaker & Mixer

Model:

Reax top /control	541-1XXXX-XX
Reax 20/4...12	541-200XX-XX
Reax 2	541-21001-XX
Duomax 1030	543-32XXX-XX
Unimax X010	54X-1XXX0-XX
Promax X020	54X-2XXXX-XX
Polymax X040	54X-4XXXX-XX
Titramax 1XXX	544-1XXXX-XX
Vibramax 1XX	544-X1200-XX
Rotamax 120	544-41200-XX
Multi Reax	545-10000-XX
Inkubator 1000	549-90010-XX

Description:

Shaker & Mixer with different working direction

The object of the declaration described above is in conformity with the relevant
UK Statutory Instruments (and their amendments):

2008 No. 1597	<i>The Supply of Machinery (Safety) Regulations</i>
2008 2016 No. 1091	<i>The Electromagnetic Compatibility Regulations 2016</i>
2012 No. 3032	<i>The Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations 2012</i>

and complies with the following technical standards :

EN ISO 12100:2010, EN 61326-1:2013, EN 61010-1:2010, EN IEC 63000:2018

UK Authorised Representative (for authorities only)

ProductIP (UK) Ltd.
8. Northumberland Av.
London WC2N 5BY

Signed for and on behalf of Heidolph Instruments GmbH & Co. KG
Walpersdorfer Straße 12, 91126 Schwabach / Germany

Schwabach, 01.09.2021

Wolfgang Jaenicke
Managing Director

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With this declaration, we confirm (according to current knowledge) that all sold laboratory devices by Heidolph Instruments GmbH & Co. KG fulfill the requirements of the EU directive 2011 / 65 / EU (RoHS) and the amended of directive 2015 / 863.

All devices are compatible with the requirement of the RoHS for the following materials:

Max. 0,01% of the weight in homogeneous material for cadmium and max. 0,1 % of the weight in homogeneous material for lead, mercury, hexavalent chromium, polybrominated biphenyl, polybrominated diphenyl ether, Di (2-ethylhexyl l) phthalate, butyl benzyl phthalate, dibutyl phthalate, diisobutyl phthalate.

In the case of individual assemblies, maximum concentrations maybe exceeded within the permissible exceptions of the Directive.

Schwabach, 22.07.2019

Stefan Peters
Vice President Marketing, Innovation & Technologie

Marcell Sarré
Vice President Quality Management & Technical Service

China RoHS Certification

China RoHS DECLARATION OF CONFORMITY

Heidolph Instruments GmbH & Co.KG has made reasonable efforts to ensure that hazardous materials and substances may not be used in its products.

In order to determine the concentration of hazardous substances in all homogeneous materials of the subassemblies, a "Product Conformity Assessment" (PCA) procedure was performed. As defined in GB/T 26572 the "Maximum Concentration Value" limits (MCV) apply to these restricted substances:

- Lead (Pb): 0.1%
- Mercury (Hg): 0.1%
- Cadmium (Cd): 0.01%
- Hexavalent chromium (Cr(VI)): 0.1%
- Polybrominated biphenyls (PBB): 0.1%
- Polybrominated diphenyl ether (PBDE): 0.1%

Environmental Friendly Use Period (EFUP)

EFUP defines the period in years during which the hazardous substances contained in electrical and electronic products will not leak or mutate under normal operating conditions. During normal use by the user such electrical and electronic products will not result in serious environmental pollution, cause serious bodily injury or damage to the user's assets.

The Environmental Friendly Use Period for Heidolph Instruments GmbH & Co.KG products is 25 years.

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This table is created according to SJ/T 11364-2014

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Note: Table applies to all products. Some of the components or parts listed above may not be part of the enclosed product.

- O: 表示该有毒有害物质在该部件所有均质材料中的含量均在GB/T 26572规定的限量要求以下。
- O: Indicates that the above mentioned hazardous substance contained in all homogeneous materials of the part is below the required limit as defined in GB/T 26572.
- X: 表示该有毒有害物质至少在该部件某一均质材料中的含量超出GB/T 26572规定的限量要求。
- X: Indicates that the above mentioned hazardous substance contained in at least one of the homogeneous materials of this part is above the required limit as defined in GB/T 26572.

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Apart from the disclosures in the above table, the subassemblies are not intentionally manufactured or formulated with lead (Pb), mercury (Hg), cadmium (Cd), hexavalent chromium (CrVI), polybrominated biphenyls (PBB), and polybrominated diphenyl ethers (PBDE).

Products manufactured by Heidolph Instruments GmbH & Co.KG may enter into further devices or can be used together with other appliances .

With these products and appliances in particular, Heidolph Instruments GmbH & Co.KG will not take responsibility for the EFUP of those products and appliances.

Schwabach, 06.08.2021

Wolfgang Jaenicke
Chief Executive Officer CEO

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Subject to change without notice. The printed version of this document is not regularly updated. The latest issue of this document can be found by visiting our homepage.

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MATERIAL CONTENT DECLARATION FOR Heidolph Instruments GmbH & Co. KG PRODUCTS
© Heidolph Instruments GmbH & Co. KG